

Human EXOSC7 protein (Recombinant) (His-Tag) (STJP018850)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	34.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23016
Gene Symbol	EXOSC7
UniProt ID	EXOS7_HUMAN
Immunogen Region	1-291aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSMASVTLS EAEKVYIVHG VQEDLRVDGR GCEDYRCVEV ETDVVSNTSG SARVKLGHTD ILVGKKAEMG TPKLEKPNEG YLEFFVDCSA SATPEFEGRG GDDLGTIAN TLYRIFNNKS SVDLKTLCIS PREHCWVLYV DVLLECGGN LFDIAISIAVK AALFNTRIPR VRVLEDEEGS KDIELSDDPY DCIRLSVENV PCIVTLCK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human EXOSC7 protein 3μg under reducing conditions Coomassie blue stain