


Human AMMECR1L protein (Recombinant) (His-Tag) (STJP018863)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	SDS-PAGE
	Host / Source	E.coli
PRODUCT PROPERTIES		
	Concentration	0.5 mg/mL
	Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
	Molecular Weight	36.9 kDa
	Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	83607
	Gene Symbol	AMMECR1L
	UniProt ID	AMERL_HUMAN
	Immunogen Region	1-310aa
	Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMGKRRCV PPLEPKLAAG CCGVKKPKLS GSGTHSHGNQ STTVPGSSSG PLQNHQHVDS SSGRENVSDL TLGPGNSPIT RMNPASGALS PLPRPNGTAN TTKNLVVTAE MCCYCFDVLY CHLYGFPQPR LPRFTNDPYP LFVTWKTGRD KRLRGCIETF SAMNLHSGLR EYTLTSALKD SRFPPLTREE LPKLFCSVSL LTNFEDAS
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human AMMECR1L protein 3 μ g under reducing conditions with Coomassie blue stain