

Human B9D2 protein (Recombinant) (His-Tag) (STJP018864)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 2mM DTT, 1mM EDTA
Molecular Weight	21.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80776
Gene Symbol	B9D2
UniProt ID	B9D2_HUMAN
Immunogen Region	1-175aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAEVHVI GQIIGASGFS ESSLFCKWGI HTGAAWKLLS GVREGQTQVD TPQIGDMAYW SHPIDLHFAT KGLQGWPRLH FQVWSQDSFG RCQLAGYGFC HVPSSPGTHQ LACPTWRPLG SWREQLARAF VGGGPQLLHG DTIYSGADRY RLHTAAGGTV HLEIGLLLRN FDRYGVEC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human B9D2 protein under reducing conditions stained with Coomassie blue