

Human PACT protein (Recombinant) (His-Tag) (STJP018871)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% Glycerol, 2mM DTT, 1mM EDTA
Molecular Weight	36.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8575
Gene Symbol	PRKRA
UniProt ID	PRKRA_HUMAN
Immunogen Region	1-313aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMQSRHR AEAPPLERED SGTFSLGKMI TAKPGKTPIQ VLHEYGMKTK NIPVYECERS DVQIHVPTFT FRVTVGDTIC TGEGTSKKLA KHRAAEAAIN ILMANASICF AVPDPLMPDP SKQPKNLNP IGSLOELAIH HGWRLPEYTL SQEGGPAHQR EYTTICRLES FMETGKGASK KQAKRNAAEK FLAKFSNISP ENHISLTN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human PACT protein 3µg under reducing conditions Coomassie blue stain