

Human DAP protein (Recombinant) (His-Tag) (STJP018874)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% Glycerol, 1mM DTT
Molecular Weight	13.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1611
Gene Symbol	DAP
UniProt ID	DAP1_HUMAN
Immunogen Region	1-102aa
Immunogen Sequence	MSSPPEG KLETKAGHPP AVKAGGMRIV QKHPHTGDTK EEKDKDDQEW ESPSPPKPTV FISGVIARGD KDFPPAAQV AHQKPHASMD KHPSPTQHI QQPRK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human DAP protein 3µg under reducing conditions with Coomassie blue stain