

Human GNG13 protein (Recombinant) (His-Tag) (STJP018884)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% Glycerol, 2mM DTT, 0.1mM PMSF
Molecular Weight	10 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51764
Gene Symbol	GNG13
UniProt ID	GBG13_HUMAN
Immunogen Region	1-64aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMEEDVDP QMKKEVESLK YQLAFQREMA SKTIPELLKW IEDGIPKDPF LNPDLMKNNP WVEKGKC

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human GNG13 protein under reducing conditions stained with Coomassie blue