

Human CD39L2/ENTPD6 protein (Recombinant) (His-Tag) (STJP018891)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	48.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	955
Gene Symbol	ENTPD6
UniProt ID	ENTP6_HUMAN
Immunogen Region	61-484aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSKWHRATA TQAFFSITRA APGARWGQQA HSPLGTAADG HEVFYGIMFD AGSTGTRVHV FQFTRPPRET PTLTHETFKA LKPGLSAYAD DVEKSAQGIR ELLDVAQQDI PFDFWKATPL VLKATAGLRL LPGEKAQKLL QKVKEVKAS PFLVGDDCVS IMNGTDEGVS AWITINFLTG SLKTPGGSSV GMLDLGGGST QIAFLPRV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human CD39L2/ENTPD6 protein under reducing conditions