

Human GLTPD1 protein (Recombinant) (His-Tag) (STJP018925)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 50% Glycerol, 1mM DTT
Molecular Weight	26.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80772
Gene Symbol	CPTP
UniProt ID	CPTP_HUMAN
Immunogen Region	1-214aa
Immunogen Sequence	MDDSETG FNLKVVLSF KQCLDEKEEV LLDPYIASWK GLVRFNLNLG TIFSFKSDV VSKLRIMERL RGGPQSEHYR SLQAMVAHEL SNRLVDLERR SHHPESGCRT VLRHLRALHW LQLFLEGLRT SPEDARTSAL CADSYNASLA AYHPWVVRRA VTVAFCTLPT REVFLFAMNV GPPEQAVQML GEALPFIQRV YNVSQK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human GLTPD1 protein under reducing conditions