

Human RPL35 protein (Recombinant) (His-Tag) (STJP018931)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	16.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11224
Gene Symbol	RPL35
UniProt ID	RL35_HUMAN
Immunogen Region	1-123aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMKIKAR DLRGKKKEEL LKQLDDLKVE LSQLRVAKVT GGAASKLSKI RVVRKSIARV LTVINQTQKE NLRKFYKGKK YKPLDLRPKK TRAMRRRLNK HEENLTKKKQ QRKERLYPLR KYAVKA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RPL35 protein 3µg under reducing conditions Coomassie blue stain