

e.coli ackA protein (Recombinant) (His-Tag) (STJP018938)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% Glycerol, 1mM DTT
Molecular Weight	45.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-400aa
Immunogen Sequence	MGSSHHHHHH SGLVPRGSH MGSMSKLV LNLGSSSLK FAIDAVNGE EYLSGLAECF HLPPEARIKWK MDGNKQEAAL GAGAAHSEAL NFIVNTILAQ KPELSAQLTA IGHRIVHGGE KYTSSVIDE SVIQGIKDAASFAPLHNPAAH LIGIEEALKS FPQLKDKNVA VFDTAFHQTM PEESYLYALP YNLYKEHGIR RYGAHGTSHF YVTQEAAKML NKPVEELN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant E. coli ackA protein under reducing conditions Coomassie blue stain