

Human Cytosolic beta-Glucosidase/GBA3 protein (Recombinant) (His-Tag) (STJP018942)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol, 1mM DTT
Molecular Weight	56.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57733
Gene Symbol	GBA3
UniProt ID	GBA3_HUMAN
Immunogen Region	1-469aa
Immunogen Sequence	MAFPAGF GWAAATAAYQ VEGGWADADGK GPCVWDTFTH QGGERVFKNQ TGDVACGSYT LWEEEDLKCIK QLGLTHYRFS LSWSRLLPDG TTGFINQKGI DYYNKIIDL LKNGVTPIVT LYHFDLPQTL EDQGGWLSEA IIESFDKYAQ FCFSTFGDRV KQWITINEAN VLSVMSYDLG MFPPGIPHFG TGGYQAAHNL IKAHAR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human Cytosolic beta-Glucosidase GBA3 protein under reducing conditions