

Human HYKK protein (Recombinant) (His-Tag) (STJP018947)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 40% Glycerol, 1mM DTT, 0.1mM PMSF
Molecular Weight	44.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	123688
Gene Symbol	HYKK
UniProt ID	HYKK_HUMAN
Immunogen Region	1-373aa
Immunogen Sequence	MSSGNYQ QSEALSKPTF SEEQASALVE SVFGLKVKV RPLPSYDDQN FHVYVSKTKD GPTEYVLKIS NTKASKNPDL IEVQNHIIMF LKAAGFPTAS VCHTKGDNTA SLVSVDGSE IKSylVRLLT YLPGRPiAEL PVSPQLLYEI GKLAakLDKT LQRfHHPKLS SLHRENfiWN LKNVPLLEKY LYALGQNRNR EIVEHV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human HYKK protein under reducing conditions with Coomassie blue stain