

Human YY1 protein (Recombinant) (His-Tag) (STJP018948)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol, 1mM DTT
Molecular Weight	47.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7528
Gene Symbol	YY1
UniProt ID	TYY1_HUMAN
Immunogen Region	1-414aa
Immunogen Sequence	MASGDTL YIATDGSEMP AEIVELHEIE VETIPVETIE TTVVGEEEEEE DDDDEDGGGG DHGGGGGGHGH AGHHHHHHHH HHHPPMIALQ PLVTDDPTQV HHHQEVLVQ TREEVVGDD SDGLRAEDGF EDQILIPVPA PAGGDDDYIE QTLVTAAAG KSGGGGSSSS GGGRVKKGGG KKSGKKSYLS GGAGAAGGGG ADPGNK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human YY1 protein 3µg under reducing conditions with Coomassie blue staining