

Human BRD1 protein (Recombinant) (His-Tag) (STJP018960)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	17.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23774
Gene Symbol	BRD1
UniProt ID	BRD1_HUMAN
Immunogen Region	556-688aa
Immunogen Sequence	EQVAMEL RLTPLTVLLR SVLDQLQDKD PARIFAQPVV LKEVPDYLDH IKHPMDFATM RKRLEAQGYK NLHEFEEDFD LIIDNCMKYN ARDTVFYRAA VRLRDQGGVV LRQARREVDV IGLEEASGMH LPERPA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human BRD1 protein 3 µg under reducing conditions
Coomassie blue stain