

Human NT5C3B protein (Recombinant) (His-Tag) (STJP018977)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol, 1mM DTT
Molecular Weight	36.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115024
Gene Symbol	NT5C3B
UniProt ID	5NT3B_HUMAN
Immunogen Region	1-300aa
Immunogen Sequence	MAEEVST LMKATVLMRQ PGRVQEIVGA LRKGGGDRQLQ VISDFDMTLS RFAYNGKRCP SSYNILDNSK IISEECRKEL TALLHHYYPI EIDPHRTVKE KLPHMVEWWT KAHNLLCQKQ IQKFQIAQVV RESNAMLREG YKTFFFNTLYH NNIFLFIFSA GIGDILEEII RQMKVFHPNI HIVSNYMDFN EDGFLQGFGK QLIHTY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NT5C3B protein 3ug under reducing conditions stained with Coomassie blue