

Mouse BMP-14/GDF-5 protein (Recombinant) (His-Tag) (STJP018984)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	16 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14563
Gene Symbol	Gdf5
UniProt ID	GDF5_MOUSE
Immunogen Region	376-495aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSAPLANRQ GKRPSKNLKA RCSRKALHVN FKDMGWDDWI IAPLEYEAFH CEGLCFPLR SHLEPTNHAV IQTLMNSMDP ESTPPTCCVP TRLSPI SILF IDSANNVVK QYEDMVVESC GCR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant mouse BMP-14 GDF-5 protein under reducing conditions