

Human DPPA5/ESG1 protein (Recombinant) (His-Tag) (STJP018991)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Phosphate buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	15.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	340168
Gene Symbol	DPPA5
UniProt ID	DPPA5_HUMAN
Immunogen Region	1-116aa
Immunogen Sequence	MGTLPAR RHIPPWVKVP EDLKDPEVFAQ VQTRLLKAIF GPDGSRIPYIEQVSKAMLEL KALESSDLTE VVYGSYLYK LRTKWMLQSM AEWHRQRQER GMLKLAEAMN ALELGPWMK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DPPA5/ESG1 protein under reducing conditions