

Human PLA2G7/PAF_AH/Lp_PLA2 protein (Recombinant) (His-Tag) (STJP019000)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	47.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7941
Gene Symbol	PLA2G7
UniProt ID	PAFA_HUMAN
Immunogen Region	47-441aa
Immunogen Sequence	AAASFGQ TKIPRGNGPY SVGCTDLMFD HTNKGTFRLR YYPSQDNDRL DTLWIPNKEY FWGLSKFLGT HWLMGNILRL LFGSMTPPAN WNSPLRPGEK YPLVVFHSHGL GAFRTLYSAL GIDLASHGFI VAAVEHRDRS ASATYYFKDQ SAAEIGDKSW LYLR TLKQEE ETHIRNEQVR QRAKECSQAL SLILDIDHGK PVKNAL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human PLA2G7 PAF-AH Lp-PLA2 protein under reducing conditions