


Human CREB3L2 protein (Recombinant) (His-Tag) (STJP019010)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli
PRODUCT PROPERTIES	
Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 50% Glycerol, 1mM DTT, 0.2M NaCl
Molecular Weight	44 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	64764
Gene Symbol	CREB3L2
UniProt ID	CR3L2_HUMAN
Immunogen Region	1-378aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMEVLESG EQGVLQWDRK LSELSEPGDG EALMYHTHFS ELLDEFNQNV LGQLNDPFL SEKSVSMEVE PSPTSPAPLI QAEHSYSLCE EPRAQSPFTH ITSDFSNDDE VESEKWYLSL DFPSTSIKTE PITDEPPPGL VPSVTLTITA ISTPLEKEEP PLEMNTGVDS SCQTIIPKIK LEPHEVDQFL NFSPKEAPVD HLHLPPTP
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CREB3L2 protein under reducing conditions