

Human TRIB3 protein (Recombinant) (His-Tag) (STJP019012)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	42.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57761
Gene Symbol	TRIB3
UniProt ID	TRIB3_HUMAN
Immunogen Region	1-358aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMRATPLA APAGSLSRKK RLEDDNLDL ERPVQKRARS GPQPRLPCL LPLSPPTAPD RATAVATASR LGPYVLEPE EGGRAYQALH CPTGTEYTCK VYPVQEALAV LEPYARLPPH KHVARPTEVL AGTQLLYAFF TRTHGDMHSL VRSRHRPEP EAAVLFRQMA TALAHCHQHG LVLRDLKLCR FVFADREKK LVLENLE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human TRIB3 protein under reducing conditions with Coomassie blue stain