

Human Epiregulin protein (Recombinant) (His-Tag) (STJP019013)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	7.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2069
Gene Symbol	EREG
UniProt ID	EREG_HUMAN
Immunogen Region	63-108aa
Immunogen Sequence	MGSSHHHHHH SGLVPRGSH MGSVSITKCS SDMNGYCLHG QCIYLVDM SQ NYCRCEVGYT GVRCEHFFL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Epiregulin protein under reducing conditions with Coomassie blue stain