

Human NMT1 protein (Recombinant) (His-Tag) (STJP019022)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4)
Molecular Weight	59.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4836
Gene Symbol	NMT1
UniProt ID	NMT1_HUMAN
Immunogen Region	1-496aa
Immunogen Sequence	MADESET AVKPPAPPLP QMMEGNGNGH EHCSDCENEE DNSYNRGGLS PANDTGAKKK KKKQKKKKEK GSETDSAQDQ PVKMNSLP AE RIQEIQAIE LFSVGQGP AK TMEEASKRSY QFWDTPVPK LGEVVNTHGP VEPDKDNIRQ EPYTLPGQFT WDALDLGDRG VLKELYTLN ENYVEDDDNM FRFDYSPEFL LWALRP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human NMT1 protein under reducing conditions stained with Coomassie blue