

Human Endophilin A3/SH3GL3 protein (Recombinant) (His-Tag) (STJP019023)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	41.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6457
Gene Symbol	SH3GL3
UniProt ID	SH3G3_HUMAN
Immunogen Region	1-347aa
Immunogen Sequence	MSVAGLK KQFHKASQLF SEKISGAEGT KLDDEFDME RKIDVTNKVV AEILSKTTEY LQPNPAYRAK LGMLNTVSKI RGQVKTTGYP QTEGLLGDCM LKYGKELGED STFGNALIEV GESMKLMAEV KDSLDIRVVKQ TFIDPLQLLQ DKDLKEIGHH LKLEGRRLD YDYKKKRVGK IPDEEVRQAV EKFEESKELA ERSMFN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Endophilin A3 SH3GL3 protein under reducing conditions