

Human Wnt-3a protein (Recombinant) (No-Tag) (STJP019029)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	37.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89780
Gene Symbol	WNT3A
UniProt ID	WNT3A_HUMAN
Immunogen Region	19-352aa
Immunogen Sequence	MSYPIWWSLA VGPQYSSLGS QPILCASIPG LVPKQLRFCR NYVEIMPSVA EGIKIGIQEC QHQFRGRRWN CTTVHDSLAI FGPVLDKATR ESAFVHAIAS AGVAFVTRTS CAEGTAAICG CSSRHQGGSPG KGWKWGGCSE DIEFGGMVSR EFADARENRP DARSAMNRHN NEAGRQAIAS HMHLKCKCHG LSGSCEVKTG WWSQPDFRAI GDFLKDKYDS ASEMVVEK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Wnt-3a protein 3 µg under reducing conditions Coomassie stain