

Human CD39L3/ENTPD3 protein (Recombinant) (His-Tag) (STJP019030)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	52.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	956
Gene Symbol	ENTPD3
UniProt ID	ENTP3_HUMAN
Immunogen Region	44-485aa
Immunogen Sequence	QIHKQEV LPPGLKYGIV LDAGSSRTTV YVYQWPAAKE NNTGVVSQTF KCSVKGSGIS SYGNNPQDVP RAFEECMQKV KGQVP SHLHG STPIHLGATA GMRLRLQNE TAANEVLESI QSYFKSQPFD FRGAQII SGQ EEGVYGWITA NYLMGNFLEK NLWHMWVHPH GVETTGALDL GGASTQISFV AGEKMDLNTS DIMQVS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CD39L3
ENTPD3 protein under reducing conditions
Coomassie stain