

Human IL-12R beta/IL12RB1 protein (Recombinant) (His-Tag) (STJP019035)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	62.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3594
Gene Symbol	IL12RB1
UniProt ID	I12R1_HUMAN
Immunogen Region	1-545aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSMEPLVTW VVPLLFLFL SRQGAACRTS ECCFQDPPYP DADSGSASGP RDLRCYRISS DRYECVSWQYE GPTAGVSHFL RCCLSSGRCC YFAAGSATRL QFSDQAGVSV LYTVTLWVES WARNTQTEKSP EVTLQLYNSV KYEPPLGDIK VSKLAGQLRM EWETPDNQVG AEVQFRHRTSP SSPWKLGDGCG PQDDDTESCL CPLEMNVA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human IL-12R beta IL12RB1 protein under reducing conditions