

Human ATP6V1F protein (Recombinant) (His-Tag) (STJP019037)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing, 50% Glycerol, 1mM DTT
Molecular Weight	15.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9296
Gene Symbol	ATP6V1F
UniProt ID	VATE_HUMAN
Immunogen Region	1-119aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAGRGKL IAVIGDEDTV TGFLGGIGE LKNRHPNFL VVEKDTTINE IEDTFRQFLN RDDIGIILIN QYIAEMVRHA LDAHQQSIPA VLEIPSKEHP YDAAKDSILR RARGMFTAED LR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ATP6V1F protein under reducing conditions stained with Coomassie blue