

Human LKB1/STK11 protein (Recombinant) (His-Tag) (STJP019042)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	51kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6794
Gene Symbol	STK11
UniProt ID	STK11_HUMAN
Immunogen Region	1-433aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSMEVDPQ QLGMFTEGEL MSVGMDTFIH RIDSTEVYQ PRRKRAKLIG KYLMGDLLGE GSYGKVKEVL DSETLCRRAY KILKKKKLRR IPNGEANVKK EIQLLRRLRH KNVIQLVDVL YNEEKQKMYM VMEYCVCGMQ EMLDSVPEKR FPVCQAHGYF CQLIDGLEYL HSQGIVHKDI KPGNLLTTG GTLKISDLGV AEALHPFA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human LKB1 STK11 protein under reducing conditions