

E.coli tnaA protein (Recombinant) (His-Tag) (STJP019053)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol, 1mM DTT
Molecular Weight	55.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-471aa
Immunogen Sequence	MENFKHL PEPFRIRVIE PVKRTTRAYR EEAIKSGMN PFLDSEDVF IDLLTDSGTG AVTQSMQAAM MRGDEAYSGS RSYIALAESV KNIFGYQYTI PTHQGRGAEQ IYIPVLIKKR EQEKGLDRSK MVAFSNYFFD TTQGHSEQING CTVRNVYIKE AFDTGVRYDF KGNFDLEGLE RGIEEVGPNN VPIYIVATITS NSAGGQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant E.coli tnaA protein under reducing conditions Coomassie blue stain