

Human LRP2BP protein (Recombinant) (His-Tag) (STJP019058)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol, 1mM DTT
Molecular Weight	42.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55805
Gene Symbol	LRP2BP
UniProt ID	LR2BP_HUMAN
Immunogen Region	1-347aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMKLTSEK LPKNPFYASV SQYAAKNQKF FQWKKEKTDY THANLVDKAL QLLKERILKG DTLAYFLRGQ LYFEEGWYEE ALEQFEEIKE KDHQATYQLG VMYYDGLGTT LDAEKGVDYM KKILDSPCPK ARHLKFAAAY NLGRAYYEGK GVKRSNEEAE RLWLIAADNG NPKASVKAQS MLGLYYSTKE PKELEKAFYW HSEACGNG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human LRP2BP protein 3µg under reducing conditions
Coomassie blue stain