

Human IP6K2 protein (Recombinant) (His-Tag) (STJP019062)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH8.0) containing 10% Glycerol
Molecular Weight	51.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51447
Gene Symbol	IP6K2
UniProt ID	IP6K2_HUMAN
Immunogen Region	1-426aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MSPAFRAMDV EPRAKGVLE PFVHQVGGHS CVLRFNETTL CKPLVPREHQ FYETLPAEMR KFTPQYKGVV SVRFEEDEDR NLCLAIYPLK GDHGIVDIVD NSDCEPKSKL LRWTTNKKHH VLETEKTPKD WVRQHRKEEK MKSHKLEEEF EWLKKSEVLY YTVEKKGNIS SQLKHYPWS MKCHQQQLQR MKENAKHRNQ YKFILLEN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human IP6K2 protein under reducing conditions Coomassie blue stain