

Human IMPAD1 protein (Recombinant) (His-Tag) (STJP019073)

GENERAL INFORMATION

Product Type	Proteins
Applications	Enzyme Activity/SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4)
Molecular Weight	37.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54928
Gene Symbol	BPNT2
UniProt ID	IMPA3_HUMAN
Immunogen Region	34-359aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSGRFSLFG LGGEPGGGAA GPAAADGGT VDLREMLAVS VLA AVRGGDE VRRVRESNVL HEKSKGKTRE GAEDKMTSGD VLSNRKMFYL LKTAFPSVQI NTEEHVDAAD QEVILWDHKI PEDILKEVTT PKEVPAESVT VWIDPLDATQ EYTEDLRKYV TTMVCVAVNG KPMLGVIHKP FSEYTAWAMV DGGSNVKARS SYNEKTPR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human IMPAD1 protein under reducing conditions stained with Coomassie blue