

Human GNG12 protein (Recombinant) (His-Tag) (STJP019088)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	10.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55970
Gene Symbol	GNG12
UniProt ID	GBG12_HUMAN
Immunogen Region	1-69aa
Immunogen Sequence	MGSSHHHHHH SGLVPRGSH MGSMSSTAS TNNIAQARRT VQQLRLEASI ERIKVS KASA DLMSYCEEHA RSDPLLIGIP TSENPFKDKK TC

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human GNG12 protein under reducing conditions stained with Coomassie blue