

Human PDGF-B protein (Recombinant) (His-Tag) (STJP019094)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	14.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5155
Gene Symbol	PDGFB
UniProt ID	PDGFB_HUMAN
Immunogen Region	82-190aa
Immunogen Sequence	CKTRTEVFE ISRRLIDRTN ANFLWPPCV EVQRCSGCCN NNRNVQCRPTQ VQLRPVQVRK IEIVRKKPIF KKATVTLEDH LACKCETVAA ARPVTRSPGG SQEQRAKTPQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human PDGF-B protein under reducing conditions Coomassie blue stain