

Human RNMTL1 protein (Recombinant) (His-Tag) (STJP019123)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	49.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55178
Gene Symbol	MRM3
UniProt ID	MRM3_HUMAN
Immunogen Region	1-420aa
Immunogen Sequence	MGSSHHHHH SSSLVPRGSH MGSMAALVRP ARFVVRPLLQ VVQAWDLAR RWVRALRRSP VKVVFPSGEV VEQKRAPGKQ PRKAPSEASA QEQRKQPLE ESASRAPSTW EESGLRYDKA YPGDRRLSSV MTIVKSRPFR EKQGKILLEG RRLISDALK A GAVPKMFFFS RLEYLKELPV DKLKGVSLIK VKFEDIKDWS DLVTPQGIMG IFAKPDHVKM TYPKTQLQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing 3µg recombinant human RNMTL1 protein under reducing conditions
Coomassie blue stain