

Human AMT protein (Recombinant) (His-Tag) (STJP019136)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 30% Glycerol, 1mM DTT
Molecular Weight	43.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	275
Gene Symbol	AMT
UniProt ID	GCST_HUMAN
Immunogen Region	29-403aa
Immunogen Sequence	AQEVLRRL TPLYDFHLAH GGKMVAFAGW SLPVQYRDSH TDShLHTRQH CSLFDVSHML QTKILGSDRV KLMESLVVGD IAE LRPNQGT LSLFTNEAGG ILDDLIVTNT SEGHLVYVSN AGCWEKDLAL MQDKVRELQN QGRDVGLEVL DNALLALQGP TAAQVLQAGV ADDLRKLPM TSAVMEVFGV SGCRVTRCGY TGEDGV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human AMT protein 3µg under reducing conditions Coomassie blue stain