

Human GID8 protein (Recombinant) (His-Tag) (STJP019154)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol, 1mM DTT
Molecular Weight	29.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54994
Gene Symbol	GID8
UniProt ID	GID8_HUMAN
Immunogen Region	1-228aa
Immunogen Sequence	MSYAEKP DEITKDEWME KLNNLHVQRA DMNRLIMNYL VTEGFKEAAE KFRMESGIEP SVDLETLDER IKIREMILKG QIQEAIALIN SLHPELLDTN RYLYFHLQQQ HLIELIRQRE TEAALEFAQT QLAEQGGEESR ECLTEMERTL ALLAFDSPEE SPFGDLLHTM QRQKVWSEVN QAVLDYENRE STPKLAKLLK LLLWAQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human GID8 protein under reducing conditions with Coomassie blue stain