

Mouse PGP9.5/Uchl1 protein (Recombinant) (His-Tag) (STJP019165)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol, 1mM DTT
Molecular Weight	27.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22223
Gene Symbol	Uchl1
UniProt ID	Uchl1_MOUSE
Immunogen Region	1-223aa
Immunogen Sequence	MQLKPME INPEMLNKVL AKLGVAGQWR FADVLGLEEE TLGSVPSPAC ALLLLFPLTA QHENFRKKQI EELKGQEVSP KVFYFMKQTIG NSCGTIGLIH AVANNQDKLE FEDGSLVKQF LSETEKLSP E DRAKCFEKNE AIQAAHDSVA QEGQCRVDDK VNFHFILFNN VDGHLIELDG RMPFPVNHGA SSEDSSLQDA AKVCRE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant mouse PGP9.5 Uchl1 protein under reducing conditions