

Human CHRNA3 protein (Recombinant) (His-Tag) (STJP019180)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	27 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1136
Gene Symbol	CHRNA3
UniProt ID	ACHA3_HUMAN
Immunogen Region	32-240aa
Immunogen Sequence	SEAEHRL FERLFEDYNE IIRPVANVSD PVIHFEVSM SQLVKVDEVN QIMETNLWLK QIWN DYKLKWNPSDYGGAEF MRVPAQKIWK PDIVLYNNAV GDFQVDDKTK ALLKYTG EVT WIPPAIFKSS CKIDVTYFPFDYQNCTMKFG SWSYDKAKID LVLIGSSMNL KDYWESGEWA IIKAPGYKHD IKYNCCEEIY PDITYS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CHRNA3 protein 3µg under reducing conditions stained with Coomassie blue