

Human Aminopeptidase P1/XPNPEP1 protein (Recombinant) (His-Tag) (STJP019181)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol
Molecular Weight	73.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7511
Gene Symbol	XPNPEP1
UniProt ID	XPP1_HUMAN
Immunogen Region	1-623aa
Immunogen Sequence	MPPKVTSE LLRQLRQAMR NSEYVTEPIQ AYIIPSGDAH QSEYIAPCDC RRAVSGFDG SAGTAIITEE HAAMWTDGRY FLQAAKQMDN NWTLMKMKGLK DTPTQEDWLV SVLPEGSRVG VDPLIIPDTY WKKMAKVLRS AGHHLIPVKE NLVDKIWTDR PERPCKPLLT LGLDYGISW KDKVADLRK MAERNV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human XPNPEP1 protein under reducing conditions stained with Coomassie blue