

Human UBE2G2 protein (Recombinant) (His-Tag) (STJP019187)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol, 1mM DTT
Molecular Weight	21 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7327
Gene Symbol	UBE2G2
UniProt ID	UB2G2_HUMAN
Immunogen Region	1-165aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMAGTALK RLMAEYKQLT LNPPEGIVAG PMNEENFFEW EALIMGPEdT CFEEGVFPAl LSFPLDYPLS PPKMRFTCEM FHPNIYPDGR VCISILHAPG DDPMGYESSA ERWSPVQSV E KILLSVVSML AEPNDESGAN VDASKMWRDD REQFYKIAKQ IVQKSLGL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human UBE2G2 protein under reducing conditions stained with Coomassie blue