

Human UBE2Q2 protein (Recombinant) (His-Tag) (STJP019189)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 30% Glycerol, 1mM DTT
Molecular Weight	45.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	92912
Gene Symbol	UBE2Q2
UniProt ID	UB2Q2_HUMAN
Immunogen Region	1-375aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMVSGLK AELKFLASIF DKNHERFRIV SWKLDLHCQ FLVPQQGSPH SLPPPLTLHC NITESYPSSS PIWFDSEDP NLTSVLERLE DTKNNNLLRQ QLKWLICELC SLYNLPKHLD VEMLDQPLPT GQNGTTEEV TSEEEEEEEEM AEDIEDLDHY EMKEEPISG KKSEDEGIEK ENLAILEKIR KTQRQDHLNG AVSGSVQA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human UBE2Q2 protein under reducing conditions with Coomassie blue stain