

Human UBCH9/UBE2E3 protein (Recombinant) (His-Tag) (STJP019196)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol, 1mM DTT
Molecular Weight	25.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10477
Gene Symbol	UBE2E3
UniProt ID	UB2E3_HUMAN
Immunogen Region	1-207aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMSDDRQR SDDSPSTSS GSSDADQRDP AAPEPEEQEE RKPSATQQKK NTKLSSKTTA KLSTSAKRIQ KELAEITLDP PPNCAGPKG DNIYEWRSI LGPPGSVYEG GVFFLDITFS SDYPPKPPKV TFRTRIYHCN INSQGVICLD ILKDNWSPAL TISKVLLSIC SLLTDCNPAD PLVGSIAQYQ LTNRAEHDRI ARQWTKRY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human UBCH9/UBE2E3 protein under reducing conditions