

## Mouse Glutathione S-Transferase Pi 1/GSTP1 protein (Recombinant) (His-Tag) (STJP019207)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	26 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">14870</a>
Gene Symbol	<a href="#">Gstp1</a>
UniProt ID	<a href="#">GSTP1_MOUSE</a>
Immunogen Region	1-210aa
Immunogen Sequence	MPPYTIV YFPVGRCEA MRMLLADQQG SWKEEVVTID TWMQGLLKPT CLYGQLPKFE DGDLTLYQSN AILRHLLGRSL GLYGKNQREA AQMDMVNDGV EDLRGKYVTL IYTNYENGKN DYVKALPGHL KPFETLLSQN QGGKAFIVGD QISFADYNLL DLLLIHQVLA PGCLDNFPLL SAYVARLSAR PKIKAFLLSSP EHVNRP

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant mouse GSTP1 protein  
3µg reducing conditions Coomassie stain