

Human Angiopoietin-like 3/ANGPTL3 protein (Recombinant) (His-Tag) (STJP019209)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	27.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27329
Gene Symbol	ANGPTL3
UniProt ID	ANGL3_HUMAN
Immunogen Region	243-460aa
Immunogen Sequence	MGSSHHHHHH SGLVPRGSH MPAECTTIYN RGEHTSGMYA IRPSNSQVVFH VYCDVISGSP WTLIQHRIDG SQNFNETWEN YKYGFGRLDG EFWLGLLEKIY SIVKQSNYVL RIELEDWKDN KHYIEYSFYL GNHETNYTLH LVAITGNVNP AIPENKDLVF STWDHKAKGH FNCPEGYSGG WWWHDEC GEN NLNGKYNKPR AKSKPERRRG LSWKSQNGRL YSIKSTKM

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human Angiopoietin-like 3 ANGPTL3 protein under reducing conditions