

Human NMNAT-1 protein (Recombinant) (His-Tag) (STJP019231)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT 1mM EDTA
Molecular Weight	36.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64802
Gene Symbol	NMNAT1
UniProt ID	NMNA1_HUMAN
Immunogen Region	1-279aa
Immunogen Sequence	MENS EKTEVLLAC GSFNPITNMH LRLFELAKDY MNGTGRTYTVV KGIISPVGDA YKKKGLIPAY HRVIMAEELAT KNSKWVEVDT WESLQKEWKE TLKVLRRHQE KLEASDCDHQ QNSPTLERPG RKRKWTETQD SSQKKSLEPK TKAVPKVKLL CGADLLESFA VPNLWKSEDI TQIVANYGLI CVTRAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NMNAT-1 protein under reducing conditions Coomassie blue stain