

Human EGFR protein (Recombinant) (hIgG-His-Tag) (STJP019271)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	CHO Cells

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	95.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1956
Gene Symbol	EGFR
UniProt ID	EGFR_HUMAN
Immunogen Region	25-645aa
Immunogen Sequence	LEEKKVCQGT SNKLTQLGTF EDHFLSLQRM FNNCEVVLGN LEITYVQRNY DLSFLKTIQE VAGYVLIALN TVERIPLLENL QIIRGNMYYE NSYALAVLSN YDANKTGLKE LPMRNLQEIL HGAVRFSNNP ALCNVESIQQW RDIVSSDFLS NMSMDFQNHG GSCQKCDPSC PNGSCWGAGE ENCQKLTII CAQQCSGRCR GKSPSDCCHN QCAAGCTGPR ESDCLVCR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human EGFR protein 3µg under reducing conditions stained with Coomassie blue