

Human S100A8 protein (Recombinant) (His-Myc-Tag) (STJP019276)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.1M NaCl
Molecular Weight	14.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6279
Gene Symbol	S100A8
UniProt ID	S10A8_HUMAN
Immunogen Region	1-93aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSFMLTELE KALNSIIDVY HKYSLIKGNF HAVYRDDLLK LLETCPQYI RKKGADVWFK ELDINTDGAV NFQEFLILVI KMGVAAHKKS HEESHKEEQK LISEEDL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human S100A8 protein 3µg reducing conditions Coomassie stain