


Human Monoglyceride Lipase/MGLL protein (Recombinant) (His-Tag) (STJP019291)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	SDS-PAGE/Enzyme Activity
	Host / Source	E.coli
PRODUCT PROPERTIES		
	Concentration	0.5 mg/mL
	Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
	Molecular Weight	36.4 kDa
	Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	11343
	Gene Symbol	MGLL
	UniProt ID	MGLL_HUMAN
	Immunogen Region	1-313aa
	Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH METGPEDPSS MPEESSPRRT PQSIPYQDLP HLVNADGQYL FCRYWKPTGT PKALIFVSHG AGEHSGRYEE LARMLMGLDL LVFAHDHVGH GQSEGERMVV SDFHVFVRDV LQHVDSMQKD YPGLPVFLLG HSMGGAAIL TAAERPGHFA GMVLISPLVL ANPESATTFK VLAAKVLNLV LPNLSLGPID SSVLSRNKTE VDIYNNDP
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Monoglyceride Lipase MGLL protein under reducing conditions