

canine IFN-gamma/IFNG protein (Recombinant) (No-Tag) (STJP019296)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	16.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	24-166aa
Immunogen Sequence	MQAMFFKEIE NLKEYFNASN PDVSDGGSLF VDILKKWREE SDKTIIQSQI VSFYKLFND FKNQIIQRS MDTIKEDMLG KFLNSSTSKR EDFLKLQIP VNDLQVQRKA INELIKVMND LSPRSNLRKR KRSQNLFRGR RASK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant canine IFN-gamma IFNG protein under reducing conditions